

CLEANING, DISINFECTING AND STERILIZING
NETTOYAGE, DÉSINFECTION ET STÉRILISATION
REINIGEN, DESINFIZIEREN UND STERILISIEREN
LIMPIEZA, DESINFECCIÓN Y ESTERILIZACIÓN
PULIZIA, DISINFEZIONE E STERILIZZAZIONE
RENGÖRING, DESINFICERING OCH STERILISERING

Reprocessing Instructions

DENTAL

ADDITIONAL SAFETY ADVICE FOR USA:

Caution! Federal (USA) law restricts this unit to sale by or on the order of a dentist

EMS 
ELECTRO MEDICAL SYSTEMS

PLEASE READ THIS BEFORE BEGINNING WORK!

Please read these operating instructions carefully as they explain all the most important details and procedures. Please pay special attention to the safety precautions.

Always keep this instruction close at hand.

To prevent injury to people and damage to property, please heed the corresponding directives. They are marked as follows:

Danger

Risk of injury

Caution

Risk of property damage or environmental harm

Please note

Useful additional information and hints

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Tillåtet

About this manual

The instructions in this manual are only applicable to the equipment it was delivered with.

Note that the English version of this manual is the master from which translations derive. In case of any discrepancy, the binding version is the English text.

We would be pleased to answer your questions or listen to your suggestions. We do, of course, provide support in case of technical problems. Please call directly your approved EMS repair center or your dealer.

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EMS

INTRODUCTION

These reprocessing instructions provide instructions for cleaning, disinfection, sterilization and packaging of EMS reusable products intended to be reprocessed in medical facilities.

The goal of reprocessing reusable products is to reduce bioburden and to achieve sterility of those products in order to eliminate the risk of product reuse related infection. Decisions regarding cleaning, disinfecting or sterilizing EMS medical and dental instruments are based on the potential risk of infection associated with their use.

It is recommended to use steam sterilization or high level disinfection with liquid chemicals.

Remember that sterilization or high-level disinfection cannot be achieved unless the elements of the assembly are cleaned first.

If there is anything in the following instructions that is not clear, do not hesitate to contact EMS in Nyon.

If you find that the manufacturer's reprocessing instructions seem to be inadequate, please inform EMS in Nyon about those inadequacies.

We encourage to report adverse events related to device reprocessing. Report such events directly to EMS in Nyon.

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CONTENTS

Reprocessing instructions for reusable products .2	
Application..... 2	
Important 2	
Preparation 4	
Basic principles..... 4	
Post-operative treatment 6	
Steps..... 6	
Cleaning/disinfection 8	
Manual procedure..... 8	
Steps for manual cleaning 12	
Steps for manual disinfection..... 14	
Inspection and maintenance..... 18	
Packaging and sterilization 20	
Service life 26	

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REPROCESSING INSTRUCTIONS FOR REUSABLE PRODUCTS

The instructions are binding for the reprocessing of all **reusable products*** of EMS. When necessary, additional product-specific instructions are included with the product to provide additional information.

*(Here after called "products")

APPLICATION

Important

Read carefully before use the operating instructions of the EMS instrument and devices with which the product will be used.

Reusable products must be replaced after the number of operations specified by manufacturer. Disposable products may not be reused.

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PREPARATION

Basic principles

It is only possible to carry out effective sterilization after the completion of effective cleaning and disinfection. Please ensure that, as part of your responsibility for the sterility of products during use, only sufficiently validated equipment and product-specific procedures are used for cleaning/disinfection and sterilization, and that the validated parameters are adhered to during every cycle.

Please also observe the applicable legal requirements in your country as well as the hygiene regulations of the hospital or clinic. This applies especially also with regard to the additional requirements for the inactivation of prions.

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POST-OPERATIVE TREATMENT

 The post-operative treatment must be carried out immediately, nevertheless within a max. of 30 minutes at the latest after the completion of the operation. Additional information, where necessary, is provided in the respective product-specific usage instructions.

Steps

1. Rinse away any surface soiling on the product with distilled deionized water or with a cleaning agent.
2. Rinse through all lumina (e.g. irrigation and aspiration connection) at least 3 times in the normal direction of flow (no back rinsing) using a disposable syringe (min. volume 50 ml) filled with distilled/deionized water applied to the back nozzle.
3. An aldehyde-free cleaning and disinfection solution that is compatible with the products (see section Service life) may also be used as an alternative rinsing solution. In such cases, it is necessary to rinse through afterwards at least 3 times with distilled, or deionized water.

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CLEANING/DISINFECTION

Manual procedure

 Subsequent cleaning/disinfection must be commenced within 2 hours at the maximum.

When selecting the cleaning and/or disinfection agent to be used, ensure that:

- These are fundamentally suitable for the cleaning and/or disinfection of the products and compatible with one another (because of the heavy soiling which may occur, the use of combined cleaning/disinfection agents are not recommended),
- The disinfection agent holds a tested effectiveness (e.g. DGHM, FDA approval or CE marking),
- The chemicals used are compatible with the products (see section Service life).

 It is absolutely essential that the concentrations and contact times specified by the manufacturer of the cleaning and disinfection agent are adhered to. Only fresh prepared solutions may be used. The disinfectant solution is not permitted to foam.

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Only sterilized or low microbe count distilled/de-ionized water (< 10 cfu/ml) can be used for all rinsing steps. Also ensure the presence of a sufficiently low endotoxin and particle concentration (e.g. Aqua purificata, as per the specifications of Pharm. Eur. or USP).

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Steps for manual cleaning

Handpiece cleaning

1. Remove any externally attached stains by brushing carefully with a soft brush or a soft cloth.
2. Rinse through all product lumina (e.g. irrigation and aspiration connection) at least 5 times in the flow direction (no back rinsing), using a disposable syringe (min. volume 50 ml) applied to the nozzles of the product.
3. Rinse thoroughly the external body of the handpiece.

Instruments and other products cleaning

1. Place the products in the cleaning solution for, as a minimum, the time and the concentration specified by the manufacturer of the cleaning/disinfection agent.
2. Remove any externally-attached soiling by brushing carefully with a soft brush or a soft cloth.
3. Rinse the products vigorously at least five times, each time with fresh distilled or deionized water (each product lumen with at least 50 ml of water). Repeat the cleaning process if the last rinsing does not run clear, or if stains are still visible on the product.

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Steps for manual disinfection

Handpiece disinfection

1. Rinse the products carefully at least five times, each time with fresh distilled or deionized water (each product lumen with at least 50 ml of water). Repeat the cleaning process if the last rinsing does not run clear, or if stains are still visible on the product.
2. Clean all the surface with alcohol based disinfection wipes according to the specification of the manufacturer of the cleaning/disinfection agent.
3. Dry the products with filtered compressed air (max. 3 bars).
4. Wrap the products directly following disinfection and drying (see section Packaging and sterilization). If necessary, repeat drying of the product in a clean place.

Instruments and accessories disinfection

1. Place the products in the disinfection solution for, as a minimum, the time specified by the manufacturer of the disinfection agent. Use an aldehyde-free disinfection solution for cleaning general and surgical instruments.

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2. Dry the products with filtered compressed air (max. 3 bars).

3. Wrap the products directly following disinfection and drying (see section Packaging and sterilization). If necessary, repeat drying of the product in a clean place.

Verification of the fundamental suitability of the products for effective manual cleaning and disinfection was provided by an independent, accredited testing laboratory for use of the following cleaning/disinfection agents:

For handpiece FD 350 (from Dürr GmbH) and for the instruments ID 212 (from Dürr GmbH). The procedure described above was adhered to as part of this.

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INSPECTION AND MAINTENANCE

 If visible stains are still visible on the product after cleaning/disinfection, the entire cleaning/disinfection procedure must be repeated. Product with visible damage, chip/flake loss, corrosion or bent must be taken out of use (no further use is permissible).

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PACKAGING AND STERILIZATION

(Fractional pre-vacuum* and steam gravity displacement procedure)

Please note your sterilizer may not operate at the validated cycle parameters. Any alternate sterilization cycles must be validated in accordance with ISO 14161:2000, "Sterilization of health care products - Biological indicators - Guidance for the selection, use, and interpretation of the results."

Do not exceed the maximum number of sterilization cycles.

Only cleaned and disinfected products are permitted to be sterilized. Prior to sterilization, the products need to be placed in a suitable sterilization container:

- In accordance with EN 868 and ISO 11607,
- Resistant to 138°C, with adequate steam permeability,
- Maintained on a regular basis.

If double, single-use sterilization packaging (double bag) is to be used, this must also comply with EN 868 and ISO 11607 and be suitable for steam sterilization (temperature resistant to 138°C with adequate steam permeability).

Use only the following listed steam sterilization procedures for sterilization; other sterilization procedures are not permissible:

- Fractional pre-vacuum procedure*,

* Fractional pre-vacuum procedure = steam sterilization with repetitive pre-vacuum.

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- Steam sterilizer in accordance with EN 13060 or EN 285 validated in accordance with AAMI TIR No. 12, DIN EN ISO 14161, DIN EN 556 (in future: EN ISO 17665) (valid commissioning and product-specific performance assessment: valid DQ, IQ, OQ and PQ),
- Maximum sterilization temperature 138°C (plus 3°C tolerance as per EN 554, ISO 11134 (in future: EN ISO 17665),
- Sterilization time at least 3 min. at 132/134°C (fractional pre-vacuum procedure). Sterilization at 134°C for a maximum of 20 minutes is permissible.

Verification of the fundamental suitability of the products for effective steam sterilization was provided by an independent, accredited testing laboratory for use of the fractionated vacuum procedure.

(EuroSelectomat from MMM Münchener Medizin Mechanik GmbH, Planegg, validated according to DIN EN 285). The procedure described above was adhered to as part of this.

The hot-air sterilization and radio-sterilization procedure may not be used (destruction of products).

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The manufacturer assumes no responsibility for the use of other sterilization procedures (e.g. ethylene oxide, formaldehyde and low temperature plasma sterilization). In such cases, please observe the respective valid standards (EN ISO 14937/ANSI AAMI ISO 14937 or the procedure-specific standard) and verify the suitability and effectiveness in principle of the procedure (if necessary, including investigations on sterilizing agent residue), taking into account the specific product geometry as part of the validation.

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SERVICE LIFE

The products have been designed for a large number of sterilization cycles. The materials used in their manufacture were selected accordingly. However with every renewed preparation for use, thermal and chemical stresses will result in ageing of the products.

If the number of permissible re-sterilization cycles is restricted, this will be pointed out in the product specific instructions.

The use of ultrasound baths and strong cleaning and disinfection fluids (alkaline pH > 9 or acid pH < 5) can reduce the life span of products. The manufacturer accepts no liability whatsoever in such cases.

The products may not be exposed to temperature above 138°C.

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EMS Corporation