


GENTLE DOES IT

THE NEW ORIGINAL AIR-FLOW[®] POWDER SOFT – A SMOOTH RIDE TO GREATER RECALL SUCCESS

200 g bottle



Developed by the inventor of the Original Air-Flow Method, the new Air-Flow Powder Soft is a perfect advocate when it comes to promoting regular supra-gingival prophylaxis. Why? Because it does not abrade the enamel.

MORE RECALLS

The Original Air-Flow Powder Soft is glycine-based and has the right consistency for particularly gentle applications, even involving sensitive teeth. Air-Flow Powder Soft does not wear down the tooth structure. Instead, it polishes dental surfaces and removes only what must be removed: dental plaque.

The Original Air-Flow Powder Soft is your key to successful recalls.



For more information >
welcome@ems-ch.com

DIRECT ORDER →

FAX +41 (0) 22 99 44 701

AIR-FLOW® POWDER SOFT

THE NEW AIR-POLISHING POWDER



Grain size ~65 µm (d50)

	Number	Order no	EUR*
1 BOX OF 4 BOTTLES WITH 200g EACH	<input type="text"/>	DV-071/A	96.-

*Recommended retail price plus applicable VAT

DEALER →

DENTAL PRACTICE STAMP/SIGNATURE →

EMS⁺
ELECTRO MEDICAL SYSTEMS

EMS
ELECTRO MEDICAL SYSTEMS SA
Chemin de la Vuarpillière 31
CH-1260 Nyon

Tel. +41 22 99 44 700
Fax +41 22 99 44 701
welcome@ems-ch.com
www.ems-dent.com