

# SWISS MASTER LIGHT®

[WWW.EMS-SWISSQUALITY.COM](http://WWW.EMS-SWISSQUALITY.COM)

EMS USA Corporation  
11886 Greenville Avenue #120  
Dallas, Texas 75243  
USA

Tel. +1 972 690 8382  
Fax +1 972 690 8981  
emsusa@emsdent.com

EMS Canada Inc.  
6525, Décarie boulevard  
Suite 205 Montréal  
(Québec) H3W 3E3  
CANADA

Tel. +1 514 736 5066  
Fax +1 514 736 5065  
info@ems-canada.ca

EMS OPTIDENT  
Valley Drive - ILKLEY  
West Yorkshire -  
ENGLAND LS29 8PB

Tel. +44 1943 60 50 50  
Fax +44 1943 60 44 22  
sales@optident.co.uk

EMS  
ELECTRO MEDICAL SYSTEMS S.A.  
Chemin de la Vuarpillière 31  
CH-1260 Nyon  
SWITZERLAND

Tel. +41 22 99 44 700  
Fax +41 22 99 44 701  
welcome@ems-ch.com  
www.ems-dent.com

**EMS**   
ELECTRO MEDICAL SYSTEMS

SWISS MASTER LIGHT®

ONE SHOT POWER LIGHT

# SWISS MASTER LIGHT®

HALOGEN LIGHT CURING UNIT

## → CURES ALL COMPOSITES


> New halogen technology for extra broad wave spectrum – highly effective in curing all existing composites.

## → BLEACHING


## → INDIVIDUAL HARDNESS CURVES

> Choose any of 3 preset programs or create your own in seconds.

DIFFERENT COMPOSITE MATERIALS, 2 mm THICK,  
WITH SWISS MASTER LIGHT®  
AT 3000 mW/cm² POLYMERIZED




ACTA STUDY  
(ACADEMIC CENTER FOR DENTISTRY AMSTERDAM)  
UNDER STANDARDIZED LAB CONDITIONS



## → EXTRA LARGE LIGHT SOURCE

> Light source diameter of 11 mm (!) covers most molars – one shot, enormous time savings.



SWISS MASTER LIGHT®  
> 11 mm Ø = 1 Shot

CONVENTIONAL LAMPS  
> 6-8 mm Ø = 4 Shots



**POWERFUL  
FAST  
RELIABLE  
SWISS  
QUALITY**

## → TOP PERFORMANCE

> More than 6x more power than conventional halogen lamps - 3000 mW/cm² due to innovative water cooling.

> Odorless and quiet even when operated at full capacity.