

**PIEZON<sup>®</sup>  
POWER  
PACK**

# **MINI MASTER<sup>®</sup>**

**EMS<sup>®</sup>**


**EMS<sup>+</sup>**  
ELECTRO MEDICAL SYSTEMS

MINI MASTER®  
PIEZON® POWER PACK

# MINI MASTER®

FROM THE INVENTOR OF THE RENOWNED PIEZON® METHOD

## PIEZON® METHOD →



### ABSOLUTELY LINEAR OSCILLATIONS →

- > for unequalled precision and gentle use of each instrument
- > Perfectly adapted for each treatment with the continuous feedback control
- > Ultrasonic and hydrodynamic effects as key to successful treatment



- > Enhanced cleaning efficiency of the irrigating solution through acoustic micro-streaming
- > Bactericidal effect and disintegration of biofilm due to ultrasound generated cavitation

### MULTIFUNCTIONAL →

all applications – from endodontics, perio, scaling and cavity preparation to restorative treatment

### STAND-ALONE →

integrated water supply with 350 ml and 500 ml bottles – just plug in and start working

### ERGONOMIC →

pushbuttons for individual power setting with permanent LED display, Endo and Dry mode, removable hose with liquid flow control knob, lightweight sterilizable handpiece, instruments with CombiTorque®, stainless-steel steribox, 2-step foot pedal (step 1: rinsing only / step 2: rinsing & ultrasound activation)

### PATIENT-FRIENDLY →

Painless treatment with electronic and feedback-control giving optimum instrument efficiency



Swiss Instruments A, P, PS



Steribox

### ALL-IN-ONE →

includes 1 sterilisable handpiece and 3 Swiss Instruments<sup>PM</sup> each with CombiTorque® in stainless-steel steribox, handpiece hose, 2 bottles, external power supply and 2-steps foot pedal

**MULTI-FUNCTIONAL  
STAND-ALONE  
ORIGINAL  
PIEZON®**



**[WWW.EMS-SWISSQUALITY.COM](http://WWW.EMS-SWISSQUALITY.COM)**

EMS  
ELECTRO MEDICAL SYSTEMS SA  
Chemin de la Vuarpillière 31  
CH-1260 Nyon

Tel. +41 22 99 44 700  
Fax +41 22 99 44 701  
[welcome@ems-ch.com](mailto:welcome@ems-ch.com)  
[www.ems-dent.com](http://www.ems-dent.com)