

# SWISS INSTRUMENTS<sup>PM</sup>

© EMS SA FA-193/EN Edition 02/2009



## PIEZON<sup>®</sup> IMPLANT CLEANING

→ THE INSTRUMENT PI



4x INSTRUMENT PI (ORD.NO DT-065A)



FILE HOLDER 120° (ORD.NO DS-010)

→ **EFFICIENT**

easily debrides implant and restoration surfaces

→ **GENTLE**

comfortable for the patient thanks to smooth surface coating and gentle tip motion

→ **SAFE**

instrument coating made of the high-tech material PEEK is fully autoclavable and safely cleans sensitive surfaces of implants and restorations; the stainless steel core makes the instrument breakproof

→ **COMPATIBLE**

fits all existing EMS Piezon<sup>®</sup> ultrasonic units

EMS  
ELECTRO MEDICAL SYSTEMS SA  
Ch. de la Vuarpillière 31  
CH-1260 Nyon  
Switzerland

Fax +41 22 99 44 701  
Tel. +41 22 99 44 700  
welcome@ems-ch.com  
www.ems-company.com

**EMS**   
ELECTRO MEDICAL SYSTEMS

# SWISS INSTRUMENTS<sup>PM</sup>



## PERIO DIAMONDS

→ INSTRUMENTS DPL3  
AND HPL3

→ **EFFICIENT**

the diamond coated Piezon<sup>®</sup> Perio instruments allow for in-depth root debridement in periodontal surgery using direct vision

→ **THOROUGH**

removes restoration overhangs and recontours narrow furcations (ondoplastics) to facilitate the patients oral hygiene efforts

→ **FLEXIBLE**

instrument HPL3 (70 µm grain diameter, EMS order no: DS-075A) is used for rough cleaning and ondoplastics, instrument DPL3 (15 µm grain diameter, EMS order no: DS-065A) is used for polishing

→ **COMPATIBLE**

fits all existing EMS Piezon<sup>®</sup> ultrasonic units and comes with its own CombiTorque<sup>®</sup>, a unique combination of torque wrench and instrument holder

EMS  
ELECTRO MEDICAL SYSTEMS SA  
Ch. de la Vuarpillière 31  
CH-1260 Nyon  
Switzerland

Fax +41 22 99 44 701  
Tel. +41 22 99 44 700  
welcome@ems-ch.com  
www.ems-company.com

**EMS**   
ELECTRO MEDICAL SYSTEMS