

SWISS MASTER LIGHT®

WWW.EMS-SWISSQUALITY.COM

EMS USA Corporation
11886 Greenville Avenue #120
Dallas, Texas 75243
USA

Tel. +1 972 690 8382
Fax +1 972 690 8981
emsusa@emsdent.com

EMS Canada Inc.
6525, Décarie boulevard
Suite 205 Montréal
(Québec) H3W 3E3
CANADA

Tel. +1 514 736 5066
Fax +1 514 736 5065
info@ems-canada.ca

EMS OPTIDENT
Valley Drive - ILKLEY
West Yorkshire -
ENGLAND LS29 8PB

Tel. +44 1943 60 50 50
Fax +44 1943 60 44 22
sales@optident.co.uk

EMS
ELECTRO MEDICAL SYSTEMS S.A.
Chemin de la Vuarpillière 31
CH-1260 Nyon
SWITZERLAND

Tel. +41 22 99 44 700
Fax +41 22 99 44 701
welcome@ems-ch.com
www.ems-dent.com

EMS 
ELECTRO MEDICAL SYSTEMS

SWISS MASTER LIGHT®

ONE SHOT POWER LIGHT

SWISS MASTER LIGHT®

HALOGEN LIGHT CURING UNIT

→ CURES ALL COMPOSITES


> New halogen technology for extra broad wave spectrum – highly effective in curing all existing composites.

→ BLEACHING


→ INDIVIDUAL HARDNESS CURVES

> Choose any of 3 preset programs or create your own in seconds.

DIFFERENT COMPOSITE MATERIALS, 2 mm THICK,
WITH SWISS MASTER LIGHT®
AT 3000 mW/cm² POLYMERIZED




ACTA STUDY
(ACADEMIC CENTER FOR DENTISTRY AMSTERDAM)
UNDER STANDARDIZED LAB CONDITIONS



→ EXTRA LARGE LIGHT SOURCE

> Light source diameter of 11 mm (!) covers most molars – one shot, enormous time savings.



SWISS MASTER LIGHT®
> 11 mm Ø = 1 Shot

CONVENTIONAL LAMPS
> 6-8 mm Ø = 4 Shots



**POWERFUL
FAST
RELIABLE
SWISS
QUALITY**

→ TOP PERFORMANCE

> More than 6x more power than conventional halogen lamps – 3000 mW/cm² due to innovative water cooling.

> Odorless and quiet even when operated at full capacity.